

## PEACE

# Pulsar Efficacy: an All-Comers rEgistry: 12-month results<sup>1</sup>

### Conclusions

- Pulsar stents showed favorable study results at 12 months in patients with an average lesion length of 11.2cm
- A Primary Patency (PP)\* of 79.5% and Freedom from Target Lesion Revascularization (FTLR) of 81.0%
- Results are similar to published data on similar lesion lengths, including 4EVER<sup>2</sup> with a PP of 81.4% and a FTLR of 89.3%
- No significant difference in PP between TASC A/B vs. TASC C/D lesions (p = 0.55) and diabetics vs. non-diabetics (p = 0.92)

### Study design

Prospective, multi-center all-comers registry.


### Endpoints

- PP\*\* at 6 and 12 months
- FTLR

### Patient characteristics

n = 118\*\*

Age, yrs <sup>†</sup>	71.9 ± 9.6	
Male	64	54.2%
Hypertension	109	92.4%
Dyslipidemia	82	69.5%
Current smoker	44	37.3%
Diabetes mellitus	37	31.4%
Obesity	48	40.7%
Renal insufficiency	14	12.1%
Ankle-brachial index <sup>‡</sup>	0.63 ± 0.16	
Walking capacity (m) <sup>‡</sup>	74.4 ± 50.8	
<b>Rutherford Classification (RC)<sup>§</sup></b>		
2	36	30.5%
3	54	45.8%
4	14	11.9%
5	13	11.0%



### Lesion characteristics

n = 118

Lesion length (mm) <sup>‡</sup>	111.5 ± 71.4	
Total occlusions	67	56.7%
Popliteal artery lesions	22	18.7%
<b>TASC classification</b>		
A	28	23.6%
B	29	24.8%
C	23	19.4%
D	38	32.2%

\*Defined as binary duplex ultrasound PSVR < 2.5 at the stented target lesion[s]

\*\*30 patients lost before 12-month follow-up [18 declined re-evaluation, 5 withdrew consent and 7 died]


†One patient without classification because of missing data

‡Data shown as mean ± SD


Results	6 months	12 months
PP	87.4%	79.5%
FTLR	93.2%	81.0%
Rutherford improvement > 1	84.4%	79.8%
Rutherford unchanged	13.8%	16.5%

PP sub-group data	6 months	12 months
Diabetics	86.0%	79.2%
Lesions > 100 mm	82.9%	78.0%
CTO	89.2%	78.1%
Popliteal segment	75.0%	71.4%
Renal insufficiency	82.8%	82.7%

## 12-month PP



## Rutherford change at 12 months




## PP and FTLR at 12 months in perspective

Studies	A.L.L. <sup>§</sup>	PP	FTLR	Total occlusions
4EVER (overall) <sup>2</sup>	7.1 cm	81.4%	89.3%	21%
PEACE <sup>1</sup>	11.2 cm	79.5%	81.0%	57%
DURABILITY <sup>3</sup>	9.6 cm	72.2%	79.1%	40%
ABSOLUTE <sup>4,5</sup>	10.1 cm	63.0%	-	37%
PEACE >100 mm <sup>3</sup>	-	78.0%	-	-
SUPERA <sup>6</sup>	7.8 cm	n/a	84.0%	-

<sup>§</sup>A.L.L. – Average Lesion Length

## Key outcomes



## Principal investigators

Dr. M Lichtenberg, Arnsberg, Germany and Prof. G Wittenberg, Bielefeld, Germany

1. Lichtenberg M et al. PEACE I all-comers registry: patency evaluation after implantation of the 4-French Pulsar-18 self-expanding nitinol stent in femoropopliteal lesions. J Endovasc Ther. 2014; 21(3):373-80. doi: 10.1583/13-4637R.1; 2. Bosiers M et al. 4-French-compatible endovascular material is safe and effective in the treatment of femoropopliteal occlusive disease: results of the 4-EVER trial. J Endovasc Ther. 2013; 20(6):746-56. doi: 10.1583/13-4437MR; 3. Bosiers M et al. Nitinol stent implantation in long superficial femoral artery lesions: 12-month results of the DURABILITY I study. J Endovasc Ther. 2009; 16(3):261-9. doi: 10.1583/08-2676.1; 4. Schillinger M et al. Balloon angioplasty versus implantation of nitinol stents in the superficial femoral artery. N Engl J Med. 2006 May 4;354(18):1879-88; 5. Schillinger M et al. Sustained benefit at 2 years of primary femoropopliteal stenting compared with balloon angioplasty with optional stenting. Circulation. 2007; 115(21):2745-9; 6. Garcia LA et al. SUPERA final 3-year outcomes using interwoven nitinol biomimetic supra stent. Catheter Cardiovasc Interv. 2017; 89(7):1259-1267. doi: 10.1002/ccd.27058.

Pulsar is a trademark or registered trademark of BIOTRONIK and/or its affiliates.